

A SHARP TOOL FOR CURRENT AND FUTURE NODES OF EUV LITHOGRAPHY

SHARP High-NA actinic Reticle Review Project

Markus Benk, Ryan Miyakawa,
Kenneth Goldberg, Patrick Naulleau

2019 International Workshop on EUV Lithography
Lawrence Berkeley National Lab, Berkeley, June 11




SHARP



SHARP

SEMICONDUCTOR HIGH-NA ACTINIC
RETICLE REVIEW PROJECT



Source: Synchrotron

Optics: Zoneplate lenses

4×NA: 0.25–0.625

Illuminator: Programmable

Navigation: 2- μm position accuracy

Throughput: up to 24 sites/hour

B


C

3 4
+
4 5

EU Y SMO
OI
LINES


How much
light do you
actually need?

Source
angular
spectrum




Aperture

Fourier
synthesis
illuminator



Zoneplate
lens




- Gold pattern on Si_3N_4 -membranes
- Magnetic mounting
- Kinematic positioning

2 mm

Zoneplates

Standard Zoneplates:

- 0.25 to 0.625 4xNA
- 6° to 10° CRA
- 5 azimuthal angles



Chip B

- Zernike Phase Contrast
- Differential Interference Contrast
- Stereoscopic imaging
- Cubic Phase Modulation

0.625 4xNA:

- 22-nm hp resolution on the mask
- 5.5 nm hp resolution wafer scale (for a 4x system)

Zoneplates

Standard Zoneplates:

- 0.25 to 0.625 4xNA
- 6° to 10° CRA
- 5 azimuthal angles

Chip B

- Zernike Phase Contrast
- Differential Interference Contrast
- Stereoscopic imaging
- Cubic Phase Modulation

Chip C


- Anamorphic zoneplates
- 0.55 4x/8xNA
- 6° CRA




ORFEO&SHARP BEST FOCUS SHIFT AT NA0.55

CENTRAL OF 7 BARS: SHARP EXPERIMENT VERSUS SIMULATION



VERTICAL




LIGHT FIELD






DARK FIELD



HORIZONTAL



Good confidence @NA0.55: ranges and trends matching



Zoneplates

Standard Zoneplates:

- 0.25 to 0.625 4xNA
- 6° to 10° CRA
- 5 azimuthal angles

Chip B

- Zernike Phase Contrast
- Differential Interference Contrast
- Stereoscopic imaging
- Cubic Phase Modulation

Chip C

- Anamorphic zoneplates
- 0.55 4x/8xNA
- 6° CRA



Zoneplates

Standard Zoneplates:

- 0.25 to 0.625 4xNA
- 6° to 10° CRA
- 5 azimuthal angles

Chip B

- Zernike Phase Contrast
- Differential Interference Contrast
- Stereoscopic imaging
- Cubic Phase Modulation

New Chip

- Anamorphic zoneplates
- 0.55 4x/8xNA
- 5.355° CRA
- Central obscuration




Pupil fill

Parametric Fill

- Quasar
- 0.33 4xNA, 6° CRA

~ 10 cycles per exposure



YAG image, 4mm below focus

3300 Flex Pupil




Liu, Proc. of SPIE 9048 90480Q (2014)
Meiling, Proc. of SPIE 8322, 83221G (2012)




Pupil fill

- Conventional
- 0.33 4xNA, 6° CRA



Liu, SPIE 90480Q (2014)


Pupil diagram



YAG image, 4mm below focus


Pupil fill

- Crosspole
- 0.33 4xNA, 6° CRA



Liu, SPIE 90480Q (2014)


Pupil diagram



YAG image, 4mm below focus


Pupil fill

- Crosspole
- 0.33 4xNA, 6° CRA



Liu, SPIE 90480Q (2014)


Pupil diagram



Modulation of flux in pupil channels


Pupil fill

- Quasar
- 0.33 4xNA, 6° CRA





Liu, SPIE 90480Q (2014)

Pupil diagram



Modulation of flux in pupil channels

3400 Flex Pupil




Little, 2017 International Workshop on EUV Lithography (2017)


Illuminator Grid

SHARP demonstration:


- 1624 pupil channels
- σ from 0.06 to 1.0



Gridded Illuminator




▪ Binary Source




▪ Modulated Source

SHARP Pupil Fill




- PFM 2-s loop time




- PFM 2-s loop time

SHARP Pupil Fill




- PFM 1-s loop time




- PFM 1-s loop time

SHARP Pupil Fill




- PFM 0.5-s loop time




- PFM 0.5-s loop time

Computed Sources




- PFM 2-s loop time



- Modeled source

Computed Sources



- PFM 1-s loop time




- Modeled source

How much
light do you
actually need?

Camera upgrade

Current CCD image sensor


- 5-s readout time



Camera upgrade

LBL detector group builds fast detectors for EUV/SXR.


- CMOS
- 500 frames per second
- virtually no readout time




Flux

SHARP image:



- 50k cts/pixel (full well exposure)
~10k photons / pixel @ 5 cts/photon
- 2000 photons/(s pixel) @ $T=5s$
~4 photons/(frame pixel) @ 500 Hz
~40 photons/(frame pixel) @ 50 Hz



- 5 photons / pixel




Autofocus



Autofocus


— 1 μm

- 80 photons / pixel
- ~25 Hz




Tilted focal plane


- 2- μm steps
trough focus



Tilted focal plane



1 μm



Continuous focus stack



Continuous focus stack



Center on
best focus


Continuous focus stack



↑ Arbitrary
step size

← Center on
best focus

Continuous focus stack



Continuous focus stack



↑ Correct
tilted plane

↑ Arbitrary
step size

← Center on
best focus

↕ Exposure
time


Continuous focus stack

Regular series: $\Delta z = 400 \text{ nm}$

— 1 μm


Current minimum step size:
 $\Delta z = 50 \text{ nm}$

Tilted focal plane:
 $\Delta z = 50 \text{ nm} / 15 \text{ lines (image)}$



#	1	8	28	56	70	56	28	8	1
---	---	---	----	----	----	----	----	---	---


z [nm]	-200	-150	-100	-50	0	50	100	150	200
----------	------	------	------	-----	---	----	-----	-----	-----





composite

 1 μm


regular



ZP interferometer




Source Optimization




- Pupil Channel

α

Source Optimization




- Pupil Channel α




- Image i_a

Source Optimization




- Pupil Channel


α



Source Optimization




- Pupil Channel α




- Image i_a

Source Optimization




- Pupil Channel α




- Image i_a

Source Optimization




- Pupil Channel
 α




- Image i_a

Source Optimization




- Pupil Channel α




- Image i_a

Source Optimization




- Pupil Channel
 α




- Image i_a

Source Optimization




- Pupil

Source Optimization




- Pupil




- Image $I = \hat{a} i_a$

Source Optimization




- Pupil




- Image $I = \hat{a} i_a$
 a

Source Optimization





- Pupil



- Image $I = \hat{a} i_a$
 a

Source Optimization



- Quasar



- 40-nm (1x) dense contacts

Source Optimization



- Quasar





- 40-nm (1x) dense contacts

Source Optimization

- Quasar





- Freeform Source





- 40-nm (1x) dense contacts

Source Optimization

- Quasar











- Freeform Source



- 40-nm (1x) dense contacts

Reconstruction of pupil fill from image



Source Optimization


Limitations:

- 30+ min
acquisition time
- Number of sources
- No through-focus
capability
- Image drift

Source Optimization

Limitations:

- 30+ min acquisition time
- Number of sources
- No through-focus capability
- Image drift



Source Optimization

Upgraded sensor with high frame rate:

- full source angular spectrum in 10s total exposure time
- Shot noise identical to regular SHARP data
- Through-focus capability
- Virtually no image drift



Source Optimization

Upgraded sensor with high frame rate:


- full source angular spectrum in 10s total exposure time
- Shot noise identical to regular SHARP data
- Through-focus capability
- Virtually no image drift



Summary

SHARP High-NA Actinic Reticle Review Project


- Emulation of imaging in EUV scanner
- Emulation of anamorphic imaging
- Fully programmable illuminator
- Diffraction-limited zoneplate lenses



Summary

Camera upgrade


- Fast autofocus
- Continuous-focus stack
- Source optimization



Thanks to
our users.

Thanks to
INTEL for funding
SHARP operations.

EUV infrastructure at
Berkeley is funded through
the EUREKA program.





EUREKA



euv  tech



ASTRILEUX
ACCELERATING THE FUTURE

CYMER®

An **ASML** company